QuantStudio 3 / QuantStudio 5

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QuantStudio Design & Analysis Software v1.2

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For Research Use Only. Not intended for any animal or human therapeutic or diagnostic use.

Information in this document is subject to change without notice.

APPLIED BIOSYSTEMS DISCLAIMS ALL WARRANTIES WITH RESPECT TO THIS DOCUMENT, EXPRESSED OR IMPLIED, INCLUDING BUT NOT LIMITED TO THOSE OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. TO THE FULLEST EXTENT ALLOWED BY LAW, IN NO EVENT SHALL APPLIED BIOSYSTEMS BE LIABLE, WHETHER IN CONTRACT, TORT, WARRANTY, OR UNDER ANY STATUTE OR ON ANY OTHER BASIS FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING BUT NOT LIMITED TO THE USE THEREOF, WHETHER OR NOT FORESEEABLE AND WHETHER OR NOT APPLIED BIOSYSTEMS IS ADVISED OF THE POSSIBILITY OF SUCH DAMAGES.

Limited Use Label License No. 358: Research Use Only

The purchase of this product conveys to the purchaser the limited, non-transferable right to use the purchased amount of the product only to perform internal research for the sole benefit of the purchaser. No right to resell this product or any of its components is conveyed expressly, by implication, or by estoppel. This product is for internal research purposes only and is not for use in commercial applications of any kind, including, without limitation, quality control and commercial services such as reporting the results of purchaser's activities for a fee or other form of consideration. For information on obtaining additional rights, please contact outlicensing@lifetech.com or Out Licensing, Life Technologies, 5791 Van Allen Way, Carlsbad, California 92008.

NOTICE TO PURCHASER: LIMITED LICENSE

This Instrument is licensed for use under U.S. Patent No. 7,670,832, owned by the University of Utah Research Foundation and licensed to Idaho Technology, Inc. No right is conveyed expressly by implication or by estoppel under any other patent claim.

TRADEMARKS

All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified.

TaqMan, Amperase, and AmpliTaq Gold are registered trademarks of Roche Molecular Systems, Inc.

RealTime StatMiner is a trademark of Integromics SL.

Microsoft, Microsoft Vista, and Internet Explorer are trademarks of Microsoft Corporation.

Adobe, Acrobat, and Reader are trademarks of Adobe Systems, Incorporated.

© 2016 Thermo Fisher Scientific Inc. All rights reserved. All.

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Properties Method Plate Run Results Export

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The screenshot displays the 'Method' configuration interface. At the top, there are tabs for 'Properties', 'Method', 'Plate', 'Run', 'Results', and 'Export'. The 'Method' tab is active. Below the tabs, the 'Experiment Method' section shows 'Volume' set to 50 μ L and 'Cover' set to 105.0 $^{\circ}$ C. The main area is a thermal profile graph with 'Hold Stage' and 'PCR Stage' sections. The graph shows temperature ramps and holds for Step 1 and Step 2. A '1' is circled next to the volume input, a '2' is circled on the first PCR hold stage, and a '3' is circled on the second PCR hold stage. At the bottom, there are 'Previous' and 'Next' buttons, with a '5' circled next to the 'Next' button. A legend at the bottom indicates data collection and pause options.

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The screenshot shows the 'Assign Targets and Samples' interface with the following numbered callouts:

- 1**: Points to a well in the plate layout grid.
- 2**: Points to the 'Quick Setup' button.
- 3**: Points to the 'Add' button in the 'Targets' table.
- 4**: Points to the 'Assign' button in the 'Targets' table.
- 5**: Points to the 'Sample Name' column in the 'Samples' table.
- 6**: Points to the 'Biological Group' column in the 'Biological Replicate Groups' table.
- 7**: Points to the 'View' button.
- 8**: Points to the 'Advanced Setup' button.
- 9**: Points to the 'Action' button in the top right corner.

Targets	Name	Reporter	Quantifier	Comments	Task	Quantity
<input type="checkbox"/>	IL8	FAM	RFQ-MDD		-	X
<input type="checkbox"/>	ACTB	FAM	RFQ-M33		-	X

Samples	Sample Name	Comments
<input type="checkbox"/>	Sample 1	
<input type="checkbox"/>	Sample 2	
<input type="checkbox"/>	Sample 3	

Biological Replicate Groups	Biological Group	Comments

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Select a target

① * Model: **Singleplex** * Select the target for this standard curve: **IL8** ②

Define the standard curve

③ * # of Points: **5** (5 Recommended)

④ * # of Replicates: **3** (3 Recommended)

⑤ * Starting Quantity: **1.0** (Enter the highest or lowest standard quantity for the standard curve)

⑥ * Serial Factor: **1.5** (Select a value from 1:10 to 10x)

5 Points X 3 Replicates = 15 Required Wells

Standard Curve Preview

1E0
2E-1
4E-2
6E-3
16E-3

Select and arrange wells for the standards

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1 **QC Summary** ▼

Flag Details

Flag:	Description	Frequency	Wells
PRFDROP	Passive reference signal changes near Cr	0	
DR1/MI11	Define acceptable delta Rn based on Cr range		
CQCONF	Low Cq confidence	0	
HIGHSD	High standard deviation in replicate group	0	
OUTLIERRG	Outlier in replicate group	1	C11

3

Flag: OUTLIERRG—Outlier in replicate group
Flag Detail: The Cr deviates significantly from Cr values in the associated replicate group.
Flagged Wells: C11
[View OUTLIERRG Troubleshooting Information](#)

Total Wells:	96	Processed Wells:	96	Manually Omitted Wells:	0	Targets Used:	1
Wells Set Up:	96	Flagged Wells:	5	Analysis Omitted Wells:	0	Samples Used:	2

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
The screenshot shows the 'Experiment Properties' dialog box in the software. The 'Experiment type' dropdown menu is highlighted with a red box, and 'Relative Standard Curve' is selected. Other visible settings include: Name (2016-01-13_134449), Instrument type (QuantStudio™ 5 System), Block type (96-Well 0.2-mL Block), Chemistry (TaqMan® Reagents), and Run mode (Standard).

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Properties	Method	Plate	Run	Results	Export
Assign Targets and Samples					
Quick Setup Advanced Setup					
Well Attributes					
Sample	New Sample ▼				
Target	New Target ▼				
Well Comments	Well Comments				
Plate Attributes					
1 Passive Reference	ROX ▼				
Reference Sample:	Sample 1 ▼				
Endogenous Control:	GAPDH ▼				

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The screenshot shows the 'Plate' configuration screen. At the top, there are tabs for 'Properties', 'Method', 'Plate', 'Run', 'Results', and 'Export'. The 'Plate' tab is selected. Below the tabs is the 'Assign Targets and Samples' section, which is divided into 'Quick Setup' and 'Advanced Setup' (the latter is active). Under 'Advanced Setup', there are three main sections: 'Targets', 'Samples', and 'Biological Replicate Groups'. The 'Biological Replicate Groups' section is highlighted with a red box and a circled '2'. It contains a table with columns for 'Biological Group' and 'Comments', and a list of two groups: 'Biological Group 1' and 'Biological Group 2'. Each group has a checkbox and an 'X' icon next to it.

Assign Targets and Samples					
Quick Setup					
Advanced Setup					
Targets					
	Name	Reporter	Quencher	Comments	Task
<input type="checkbox"/>	IL8	FAM	NFQ-MGB		-
<input checked="" type="checkbox"/>	GAPDH	FAM	NFQ-MGB		-
Samples					
	Sample Name	Comments			
<input type="checkbox"/>	Sample 1				X
<input type="checkbox"/>	Sample 2				X
<input type="checkbox"/>	Sample 3				X
<input checked="" type="checkbox"/>	Sample 4				X
Biological Replicate Groups					
	Biological Group	Comments			
<input type="checkbox"/>	Biological Group 1				X
<input checked="" type="checkbox"/>	Biological Group 2				X

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Cr Settings Flag Settings **1** Relative Quantification Settings Advanced Settings Standard Curve Settings

Analysis Type
Select the type of analysis to perform.
 Multiplex Singleplex

2 **Select Reference(s)**
Select the biological replicate group and/or sample to use as the reference(s) for this experiment.
Biological Replicate Group Reference: Biological Replicate Group 1 Reference Sample: Unk_800

Select Endogenous Control(s)
Select the target to use as the endogenous control for this experiment. To use multiple endogenous controls, select Multiple Controls from the drop-down list, then select the controls as necessary.
Endogenous Control(s): IPC

Outlier Rejection
Select to reject replicates with ΔCt values less than or equal to the value entered below. These analysis settings apply only to multiplex reactions.
 Reject Replicates with specified ΔCt
 $\Delta Ct \leq$ 1.0

RQ Min/Max Calculations
Select the algorithm to use to determine RQ Min and Max values (error bars).
 Confidence Level: 95.0 % Standard Deviations: 1

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Properties	Method	Plate	Run	Results	Export
Experiment Properties Save					
Name	2016-01-13_150903				
Barcode	Barcode (optional)				
User name	User name (optional)				
Instrument type	QuantStudio™ 5 System				
Block type	96-Well 0.2-ml Block				
Experiment type	Comparative Ct ($\Delta\Delta C_T$)				
Chemistry	TaqMan® Reagents				
Run mode	Standard				
Manage chemistry details					

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The screenshot shows the 'Plate' tab in a software interface. The 'Assign Targets and Samples' section is active, with 'Advanced Setup' selected. It contains two main sections: 'Targets' and 'Samples'. The 'Targets' section has a table with columns: Name, Reporter, Quenchor, Comments, Task, and Quantity. The 'Samples' section has a table with columns: Sample Name and Comments. Below these is the 'Biological Replicate Groups' section, which is circled in red and has a '2' in a circle next to it. It contains a table with columns: Biological Group and Comments. The 'Add' button is visible in the top right of the 'Biological Replicate Groups' section.

Assign Targets and Samples						
Quick Setup		Advanced Setup				
- Targets						
<input type="checkbox"/>	IL8	FAM	NFQ-MGB		-	X
<input type="checkbox"/>	GAPDH	FAM	NFQ-MGB		-	X
- Samples						
<input type="checkbox"/>	Sample 1					X
<input type="checkbox"/>	Sample 2					X
<input type="checkbox"/>	Sample 3					X
<input type="checkbox"/>	Sample 4					X
- Biological Replicate Groups						
<input type="checkbox"/>	Biological Group 1					X
<input type="checkbox"/>	Biological Group 2					X

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

Cr Settings Flag Settings **1** Relative Quantification Settings Advanced Settings

Analysis Type
Select the type of analysis to perform.
Multiplex Singleplex

2 **Select Reference(s)**
Select the biological replicate group and/or sample to use as the reference(s) for this experiment.
Biological Replicate Group Reference: Reference Sample:

Select Endogenous Control(s)
Select the target to use as the endogenous control for this experiment. To use multiple endogenous controls, select Multiple Controls from the drop-down list, then select the controls as necessary.
Endogenous Control(s):

3 **Efficiency**
Enter percentage values between 1 and 150%

	Target	Efficiency (%)
IPC		100.0
RNaseP		100.0

Outlier Rejection
Select to reject replicates with ΔCr values less than or equal to the value entered below. These analysis settings apply only to multiplex reactions.
 Reject Replicates with specified ΔCr
ΔCr ≤

RQ Min/Max Calculations
Select the algorithm to use to determine RQ Min and Max values (error bars).
 Confidence Level: % Standard Deviations:

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The screenshot shows the 'Experiment Properties' dialog box. At the top, there are tabs for 'Properties', 'Method', 'Plate', 'Run', 'Results', and 'Export'. The 'Properties' tab is active. The form contains the following fields:

- Name:** 2016-01-13_153505 (circled 1)
- Barcode:** Barcode - optional
- User name:** User name - optional
- Instrument type:** QuantStudio™ 5 System (dropdown, circled 2)
- Block type:** 96-Well 0.2-mL Block (dropdown, circled 3)
- Experiment type:** Genotyping (dropdown, circled 4)
- Chemistry:** TaqMan® Reagents (dropdown, circled 5)
- Run mode:** Standard (dropdown, circled 6)

At the bottom right, there is a 'Next' button (circled 7) and a 'Save' button.

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The screenshot shows the 'Method' configuration screen in the QuantStudio software. The interface is divided into several sections:

- Properties:** Method (selected), Plate, Run, Results, Export.
- Experiment Method:** Includes 'Action' and 'Save' buttons.
- Volume and Cover:** Volume is set to 25 μ L, and Cover is set to 105.0 °C.
- Thermal Cycle Graph:** Shows a graph with four stages:
 - Pre-Read Stage:** Step 1, 1.6 °C/s, 25.0 °C, 00:30.
 - Hold Stage:** Step 1, 95.0 °C, 10:00, 1.6 °C/s.
 - PCR Stage:** Step 1, 95.0 °C, 00:15, 1.6 °C/s; Step 2, 60.0 °C, 01:00, 1.6 °C/s.
 - Post-Read Stage:** Step 1, 25.0 °C, 00:30, 1.6 °C/s.
- Include Checkboxes:** Pre-PCR Read, Amplification, Post-PCR Read.
- Legends:** Data Collection On, Data Collection Off, Pause On, Pause Off, Advanced Settings, V VeriFlex.
- Navigation:** Previous and Next buttons, with the Next button highlighted by a circled '5'.

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
1 Call Settings Ct Settings Flag Settings Advanced Settings

2 Data Analysis Settings

Analyze Data from Post-PCR Read Only Analyze Data from Pre-PCR Read and Post-PCR Read

Analyze Real-Time Rn Data Analyze Real-Time Rn - Median (Rna to Rnb)

Analyze Real-Time dRn Data

Default Call Settings
Default cal settings are used to make allele calls for SNP assays without custom settings. To edit the default settings, click Edit Default Settings.

Autocaller Enabled: Yes Keep Manual Calls from Previous Analysis: No Quality Value: 95 **3** Edit Default Settings

SNP Assay	Analysis Type	Autocaller	Keep Manual Calls	Quality Value	Call Settings for C__29086771_20
C__29086771_20	Default	Yes	No	95	Apply Call Settings: <input checked="" type="checkbox"/> Default Settings
C__10048053_10	Default	Yes	No	95	<input checked="" type="checkbox"/> Autocaller Enabled <input type="checkbox"/> Keep Manual Calls from Previous Analysis Quality Value <input type="text" value="95.0"/>

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
1 QC Summary ▾

Flag Details

Flag:	Description	Frequency	Wells
BADROX	Bad passive reference signal	0	
OFFSCALE	Fluorescence is offscale	0	
NOSIGNAL	No signal in well		
PCFAIL	Positive control failed	0	
SMCLUSTER	Small number of samples in cluster	2	C__10048053_10, C__29086771_20
AMPNC	Amplification in negative control	0	
DRUMIN	Define acceptable data Rn based on Ct range		

2

3

Flag: SMCLUSTER—Small number of samples in cluster
Flag Detail: The number of data points in the cluster is less than the threshold.
Flag Criteria: Number of data points in the cluster <= 2.0
Wells Containing SNP Assays: A1, A2, A3, A4, A5, A6, A7, A8, A9, A10, A11, A12, B1, B2, B3, B4, B5, B6, B7, B8, B9, B10, B11, B12, C1, C2, C3, C4, C5, C6, C7, C8, C9, C10, C11, C12, D1, D2, D3, D4, D5, D6, D7, D8, D9, D10, D11, D12, E1, E2, E3, E4, E5, E6, E7, E8, E9, E10, E11, E12, F1, F2, F3, F4, F5, F6, F7, F8, F9, F10, F11, F12, G1, G2, G3, G4, G5, G6, G7, G8, G9, G10, G11, G12, H1, H2, H3, H4, H5, H6, H7, H8, H9, H10, H11, H12
Flagged SNP Assays: C__10048053_10, C__29086771_20
[View SMCLUSTER Troubleshooting Information](#)

Total Wells:	96	Processed Wells:	96	Manually Omitted Wells:	0	SNP Assays Used:	2
Wells Set Up:	96	Flagged Wells:	0	Analysis Omitted Wells:	0	Samples Used:	7

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For Research Use only. Not for use in diagnostic procedures. © 2016 Thermo Fisher Scientific Inc. All rights reserved.
TaqMan is a registered trademark of Roche Molecular Systems, Inc.
Microsoft and Windows are registered trademarks of Microsoft Corporation.